

Nigericin, Sodium Salt 尼日利亚菌素 (钠盐)

| 产品编号 | 产品名称 | 包装规格 |
|--------------|------------------------------------|------|
| NBS7623-5mg | Nigericin, Sodium Salt 尼日利亚菌素 (钠盐) | 5mg |
| NBS7623-10mg | Nigericin, Sodium Salt 尼日利亚菌素 (钠盐) | 10mg |
| NBS7623-25mg | Nigericin, Sodium Salt 尼日利亚菌素 (钠盐) | 25mg |

产品简介:

尼日利亚菌素 (Nigericin), 来源于吸水链霉菌 (*Streptomyces hygroscopicus*) 代谢产物的一种聚醚抗生素, 常常用作钾离子载体, 促进穿透生物膜的 K^+/H^+ 交换, 作用机制类似于缬氨霉素 (货号: NBS7624-10mg), 其携带单价阳离子跨膜的特异性为:

$K^+ > Rb^+ > Cs^+ > Na^+$ 。还可用作 Pb^{2+} 离子载体, 但对其他二价阳离子未表现出载运活性。

作为一种有效的研究工具, 能够激活 ATPase 酶活性, 并且破坏线粒体膜电位。通过加强离子的膜扩散能力来杀死细菌。类似于 ATP, 尼日利亚菌素通过缝隙连接蛋白 1 (Pannexin-1) 依赖性的方式, 刺激胞内 K^+ 外流, 诱导 NLRP3 炎症小体组装, 凋亡相关蛋白 Caspase-1 成熟, 进而刺激 IL-1 β 的加工和释放。

本品以尼日利亚菌素的钠盐形式提供, CAS NO. 28643-80-3, 普遍用来研究各种生物系统中膜电位变化产生的生理后果。也常用作 NLRP3 炎症小体激活剂, 研究炎症小体应对不同刺激物和危险信号的应答方式, 从而开发靶向多种免疫疾病的治疗手段。

产品特性

1) CAS NO: 28643-80-3

2) 化学名: (R)-2-((2R,3S,6R)-6-(((2S,4R,5R,7R,9R,10R)-2-((2S,2'R,3'S,5R,5'R)-5'-((2S,3S,5R,6R)-6-hydroxy-6-(hydroxymethyl)-3,5-dimethyltetrahydro-2H-pyran-2-yl)-2,3'-dimethyloctahydro-[2,2'-bifuran]-5-yl)-9-methoxy-2,4,10-trimethyl-1,6-dioxaspiro[4.5]decan-7-yl)methyl)-3-methyltetrahydro-2H-pyran-2-yl)propanoate, monosodium salt

3) 英文同义名: Antibiotic K178, Antibiotic X464, Azalomycin M, Helexin C, Polyetherin A, BRN 1696755, NSC 292567

4) 中文同义名: 抗生素 K178, 抗生素 X464, 阿扎霉素 M, 常春藤皂甙 C, 聚醚菌素 A

5) 分子式: $C_{40}H_{67}NaO_{11}$

- 6) 分子量: 746.94 g/mol
- 7) 纯度: $\geq 98\%$
- 8) 外观: 白色至类白色固体
- 9) 溶解性: 溶于乙醇 (20 mg/ml), 几乎不溶于水

保存条件:

-20°C 干燥避光保存, 3 年有效。

产品使用:【源自文献, 仅作参考】

文献 1, Sawatari Y, Yokota A. Diversity and mechanisms of alkali tolerance in lactobacilli [published correction appears in Appl Environ Microbiol. 2007 Aug;73(16):5385]. PMID: 17449704

体外研究:

细胞类型 (Cell type): Four strains of lactobacilli and *E. hirae* JCM 8729T

实验方法 (Assay): To evaluate the effect of Δ pH on alkali tolerance, growth experiments using four strains of lactobacilli and *E. hirae* JCM 8729T were conducted in the presence or absence of nigericin at pH values between 6.5 and 9.0. Nigericin in ethanol solution was added at a final concentration of 7 μ M. Control cultures received 0.1% (vol/vol) ethanol. Specific growth rates were calculated during the exponential growth phase as described above.

文献 2, Deng CC et al. Nigericin selectively targets cancer stem cells in nasopharyngeal carcinoma. Int J Biochem Cell Biol. 2013 Sep;45(9):1997-2006. PMID: 23831840

体内研究:

动物模型 (Animal Model): Nude mice injected with S19 cells

实验方法 (Assay): Mice are randomly divided into four groups with six animals each (control, DDP, Nigericin and DDP combined with Nigericin). DDP (2.5 mg/kg) is injected intraperitoneally for five continuous days and nigericin (4 mg/kg) is administrated intraperitoneally every two days.

给药途径 (Administration): Intraperitoneal (i.p.) injection

注意事项:

1. 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

其他离子载体产品:

| 产品编号 | 产品名称 | CAS NO | 包装规格 |
|--|---|--------------|-------------------|
| 1) 钙离子载体 (Calcium Ionophore/ Ca^{2+} Ionophore) | | | |
| <u>NBS7634-1mg</u> | <u>lonomycin 离子霉素 (游离酸)</u> | 56092-81-0 | 1mg |
| <u>NBS7635-1mg</u> | <u>lonomycin 离子霉素 (钙盐)</u> | 56092-82-1 | 1mg |
| <u>NBS7632-1mg</u> | <u>A23187 钙离子载体</u> | 52665-69-7 | 1mg |
| <u>NBS7633-1mg</u> | <u>4-bromo-A23187 钙离子载体</u> | 76455-82-8 | 1mg |
| 2) 钠离子载体 (Sodium Ionophore/ Na^{+} Ionophore) | | | |
| <u>NBS7621-500ug</u> | <u>SQI-Et (Na) 钠离子载体</u> | N/A | 500 μg |
| <u>NBS7622-500ug</u> | <u>SQI-Pr (Na) 钠离子载体</u> | 1022595-16-9 | 500 μg |
| <u>NBS7620-1ml</u> | <u>莫能霉素溶液 (1,000\times)</u> | 22373-78-0 | 1ml |
| 3) 钾离子载体 (Potassium Ionophore/ K^{+} Ionophore) | | | |
| <u>NBS7623-5mg</u> | <u>Nigericin 尼日利亚菌素钠盐</u> | 28643-80-3 | 5mg |
| <u>NBS7626-5mg</u> | <u>Nonactin 无活菌素</u> | 6833-84-7 | 5mg |
| <u>NBS7625-25mg</u> | <u>Potassium Ionophore II (BB15C5) 钾离子载体 II</u> | 69271-98-3 | 25mg |